



Bottle PE / 1000 mL / GL 45

6.1608.040

For Exchange Units. Bottle for auxiliary solutions.

Phụ kiện tùy chọn

Order no.	Description
6.1602.150	Bottle attachment / GL 45 - 3 x UNF 10/32

For the connection of 1/16 in. capillaries. Used in conjunction with MSM auxiliary solutions and in Metrohm Inline Sample Preparation.



6.1608.023

Amber glass bottle / 1000 mL / GL 45

For Exchange Units. Bottle for auxiliary solutions.


