



Liquid Handling Station, left

6.2841.120

Combined rinsing and dilution station for sample changers with a Swing Head. To be assembled to the left of the tower.

Phạm vi giao hàng 6.2841.120

Qt.	Order no.	Description
1 PCS	6.2621.140	Hexagon key 2.5 mm



1 PCS

6.1808.050

Tubing adapter olive / M8 outer

1 M8 outer thread and 1 tubing olive. E.g. for thermostat jacket of exchange units and stability measuring instruments.



1 PCS

6.1801.120

PVC tubing / 4 mm / 6 mm / 2 m



1 PCS

6.2621.030

Hexagon key 4 mm



2 PCS

6.2744.060

Threaded stopper

For UNF 10/32. Stopper for IC, e.g. for the sealing of columns.



Phụ kiện tùy chọn

Order no.	Description	
6.2762.100	Mixing vessel for Liquid Handling Station Replacement mixing vessel for combined rinsing and dilution stations.	
6.2762.200	Transparent cover for the 6.2762.100 mixing vessel, PMMA Transparent cover for the 6.2762.100 mixing vessel, for visual inspection in the mixing vessel.	
2.741.0010	741 Magnetic Stirrer Magnetic stirrer for sample changers. The magnetic stirrer is located beneath the sample rack.	