



Equipment for pH Adjustment at the 845 Eluent Synthesizer

6.5329.000

pH adjustment to Eluent Synthesizer. Includes pH electrode and additional accessories.

Phạm vi giao hàng 6.5329.000

Qt.	Order no.	Description
1 PCS	6.1446.000	SGJ stopper / B-14/(15)






Eluent bottle head to connect up to 4 tubings with M6 connectors. For eluent bottle of the Eluent synthesizer.



Connection eluent bottle at the eluent Synthesizer - Titrino



Additional remote interface for the connection of devices that can be controlled via remote lines. With permanently attached cable.



1 PCS

6.2726.090

Antidiffusion valve

For antidiffusion buret tips 6.1541.XXX, 6.1543.130, 6.1543.140,
6.1543.150, 6.1543.160, 6.1543.200



1 PCS

6.3026.210

Exchange Unit 10 mL

Exchange unit with integrated data chip with 10 mL glass cylinder and light protection. PCTFE/PTFE flat stopcock, FEP tubing connection, antidiffusion buret tip and standard amber glass reagent bottle.

