



## IC equipment : PCR with dosing unit

6.5330.400

The IC equipment: PCR with dosing unit makes it possible to convey a PCR reagent with the 800 Dosino, for example for the photometric determination of UV/VIS-active components after post-column derivatization.

Phạm vi giao hàng 6.5330.400

| Qt.   | Order no.  | Description                          |
|-------|------------|--------------------------------------|
| 2 PCS | 6.1608.030 | Round glass bottle / 1000 mL / GL 45 |



3 PCS

6.1805.110

FEP tubing / M6 / 80 cm

With light and kink protection



---

2 PCS

6.1619.000

Adsorber tube for 807 Dosing Unit



---

1 PCS

6.2744.010

Pressure screw 5x

With UNF 10/32 connection. For the connection of PEEK capillaries



---

1 PCS

6.2744.080

M6 thread / UNF 10/32 adapter

With M6 outer thread. For IC instruments.



2 PCS

6.1819.100

FEP aspiration tubing to 2 L bottle

Aspiration tubing to bottle attachment with M6 connector. For the production of eluents.



2 PCS

6.1602.160

Bottle attachment GL 45 for eluents and reagents

For eluent, reagent, and rinsing bottles, with connectors for the adsorber tube and the aspiration tubing.



1 PCS

6.2744.200

Adapter UNF 10/32 outer / M6 inner

Connecting 1/16 in. PEEK capillaries and transfer tubing (M6)



2 PCS

6.1624.000

Adapter SGJ 14 for 6.1619.XXX Adsorber tube

Adapter SGJ 14. Part of 6.1619.010 for Dilution Sample Processor.



1 PCS

6.1831.030

PEEK capillary 0.75 mm i.d. / 3 m

For all IC components.



---

1 PCS

6.2057.250

Dosino holder for PCR

Dosino holder for IC instruments for PCR Dosino.



---

1 PCS

6.3033.150

Dosing unit PCR 5 mL

5 mL dosing unit specially suitable for conveying solutions for post-column derivatization (PCR) in IC.

