



IC equipment: Thin-Layer cell: Carb (Au,Pd)

6.5337.210

Equipment comprising Thin-Layer cell with additional accessories, with electrodes for carbohydrate analysis with a gold working electrode and a Pd reference electrode.

Phạm vi giao hàng 6.5337.210

Qt.	Order no.	Description
1 PCS	6.2621.140	Hexagon key 2.5 mm



1 PCS

6.2744.014

Pressure screw 2x

With UNF 10/32 connection. For the connection of PEEK capillaries



1 PCS

6.1831.010

PEEK capillary 0.25 mm i.d. / 3 m

For all IC components.



1 PCS

6.2802.000

Polishing set for solid-state electrodes

Polishing set with 1 polishing cloth and about 2 g aluminum oxide powder (grain size 0.3 μm)



1 PCS

6.1257.820

50 μm spacer to Thin-Layer cell

50 μm spacer for the amperometric Thin-Layer cell, 3 pcs.



1 PCS

6.1257.210

Au working electrode, 3 mm

Working electrode made of gold for amperometric detection, 3 mm diameter.



1 PCS

6.1257.840

25 μ m spacer to Thin-Layer cell

25 μ m spacer for the amperometric Thin-Layer cell, 3 pcs.



1 PCS

6.1257.740

Pd reference electrode to 6.1257.xxx

Palladium reference electrode for amperometric detection.

