



VA electrode equipment with Multi-Mode Electrode pro for Professional VA instruments

6.5339.030

Complete electrode set for polarographic and voltammetric determinations. Includes Multi-Mode Electrode pro, reference electrode, platinum auxiliary electrode, measuring vessel, stirrer, electrolyte solution and additional accessories for setting up and operating the Multi-Mode Electrode.

Phạm vi giao hàng 6.5339.030

Qt.	Order no.	Description
1 PCS	6.0343.100	Separate Pt Rod Electrode
		Platinum auxiliary electrode for voltammetry and CVS.



1 PCS

6.0728.120

Ag/AgCl reference electrode

Inner reference system with electrolyte $c(\text{KCl}) = 3 \text{ mol/L}$. For use with electrolyte vessel made of plastic 6.1245.010.



1 PCS

6.1204.500

Stirrer for Professional VA instruments

Stirrer for Professional VA instruments, complete, metal-free.



1 PCS

6.1226.030

Mercury drop capillary (non-silanized) for MME

Non-silanized mercury drop capillary (10 pieces)



1 PCS

6.1226.050

Mercury drop capillary (silanized) for MME

Silanized mercury drop capillary (10 pieces)



1 PCS

6.1245.010

Electrolyte vessel made of plastic

Electrolyte vessel made of plastic with ceramic diaphragm for double junction reference electrode.




1 PCS

6.1246.120

Multi-Mode Electrode pro

Mercury electrode for voltammetry. Can be operated as DME, SMDE or HMDE.



1 PCS

6.1247.020

Needles for MME

Needles for Multi Mode Electrode (3 pieces).



1 PCS

6.1415.210

Titration vessel / 10-90 mL



1 PCS

6.2308.020

Electrolyte 3 mol/L KCl (250 mL)

Electrolyte solution $c(\text{KCl}) = 3 \text{ mol/L}$, (for Ag/AgCl reference systems)



1 PCS

6.2406.000

Mercury drop catcher

Silver wire (to capture mercury drops by amalgamation)



1 PCS

6.2711.030

Tray for filling the MME

For 884 Professional VA, 797 VA Computrace, 757 VA Computrace and 647, 663, 694, 747 VA Stands



1 PCS

6.2816.020

Syringe 10 mL with Luer connection

For various applications in IC and VA



Needle for filling in mercury into the Multi-Mode Electrode and for sample addition in Karl Fischer titration.

