



CVS electrode equipment with 2 mm platinum electrode for Professional VA/CVS instruments for intensive use

6.5339.060

Complete electrode set for the determination of organic additives in electroplating baths with CVS (Cyclic Voltammetric Stripping). Contains rotating platinum working electrode (electrode diameter: 2 mm), driving axle with mercury contact, maintenance-free gel reference electrode, platinum auxiliary electrode, measuring vessels and KCl electrolyte solution.

Phạm vi giao hàng 6.5339.060

Qt.	Order no.	Description
-----	-----------	-------------

1 PCS

6.2308.020

Electrolyte 3 mol/L KCl (250 mL)

Electrolyte solution $c(\text{KCl}) = 3 \text{ mol/L}$, (for Ag/AgCl reference systems)



1 PCS

6.1415.210

Titration vessel / 10-90 mL



1 PCS

6.1415.250

Titration vessel / 50-150 mL



1 PCS

6.0343.100

Separate Pt Rod Electrode

Platinum auxiliary electrode for voltammetry and CVS.



1 PCS

6.0730.100

Complete LL-Ag/AgCl(gel) reference electrode

Maintenance-free Double-Junction Ag/AgCl reference electrode for voltammetry. Made of glass with electrolyte $c(\text{KCl}) = 3 \text{ mol/L}$ (gel), for pH 0 - 8. Suitable for use in continuous operation in automated systems and online instruments.



1 PCS

6.1204.610

2 mm platinum electrode tip for CVS

Metal electrode tip made of platinum (Pt), diameter of electrode disk 2 mm $\pm 0.05 \text{ mm}$, polished surface, shaft made of glass, M3 thread connection. Concentricity error: 0.3 mm. Temperature range: 0 - 50 °C.



1 PCS

6.1204.520

Driving axle for rotating disk electrode (RDE) with mercury contact

Driving axle for rotating disk electrode (RDE) with titanium axle and mercury contact for Professional VA/CSV instruments, M3 thread.

