



CVS electrode equipment with 2 mm platinum electrode, gel reference electrode and steel auxiliary electrode for Professional VA/CVS instruments

6.5339.110

Complete electrode set for the determination of organic additives in electroplating baths with CVS (Cyclic Voltammetric Stripping). Contains rotating platinum working electrode (electrode diameter: 2 mm), driving axle with mercury contact, maintenance-free gel reference electrode, steel auxiliary electrode, measuring vessels and KCl electrolyte solution.

Phạm vi giao hàng 6.5339.110

Qt.	Order no.	Description
-----	-----------	-------------

1 PCS

6.0343.110

Separate steel electrode

Steel auxiliary electrode for CVS



1 PCS

6.0730.100

Complete LL-Ag/AgCl(gel) reference electrode

Maintenance-free Double-Junction Ag/AgCl reference electrode for voltammetry. Made of glass with electrolyte $c(\text{KCl}) = 3 \text{ mol/L}$ (gel), for pH 0 - 8. Suitable for use in continuous operation in automated systems and online instruments.



1 PCS

6.1204.520

Driving axle for rotating disk electrode (RDE) with mercury contact

Driving axle for rotating disk electrode (RDE) with titanium axle and mercury contact for Professional VA/CSV instruments, M3 thread.



1 PCS

6.1204.610

2 mm platinum electrode tip for CVS

Metal electrode tip made of platinum (Pt), diameter of electrode disk $2 \text{ mm} \pm 0.05 \text{ mm}$, polished surface, shaft made of glass, M3 thread connection. Concentricity error: 0.3 mm. Temperature range: 0 - 50 °C.



1 PCS

6.1415.210

Titration vessel / 10-90 mL



1 PCS

6.1415.250

Titration vessel / 50-150 mL



1 PCS

6.2308.020

Electrolyte 3 mol/L KCl (250 mL)

Electrolyte solution $c(\text{KCl}) = 3 \text{ mol/L}$, (for Ag/AgCl reference systems)

