



VA electrode shaft accessories with SPE electrode shaft for VA Computrace

6.5340.030

Accessory equipment for use with *screen-printed electrodes* (SPE). Contains electrode shaft for screen-printed electrodes, measuring vessel, stirrer, and additional accessories. Without electrodes.

Phạm vi giao hàng 6.5340.030

Qt.	Order no.	Description
1 PCS	6.1204.200	Stirrer for VA Stand
		Stirrer for VA Stand, complete, metal-free



1 PCS

6.1241.090

Electrode shaft for screen-printed electrodes (SPE)

Holder for *screen-printed electrodes* (SPE) with plug contacts for connecting the electrode cables for working, reference, and auxiliary electrodes. When used with screen-printed electrodes, the electrode shaft provides a complete sensor for voltammetry. Whereas the screen-printed electrodes must be replaced at regular intervals, the electrode shaft can continue to be used.




1 PCS

6.1244.020

Driving belt

For Rotating Disk Electrodes and stirrer in VA Stands (3 pieces).



1 PCS

6.1244.060

Silicone seal for SPE electrode shaft (10 pieces)

Seal made of silicone for SPE electrode shaft 6.1241.090 and measuring head (SPE) 6.01256.020, 10 pieces.



1 PCS

6.1415.150

Titration vessel / 5-70 mL



1 PCS

6.1446.000

SGJ stopper / B-14/(15)



4 PCS

6.1446.040

Threaded stopper / M6



1 PCS

6.1808.000

Connection bushing M6



1 PCS

6.1819.000

Gas inlet tube

In lid of measuring vessel in VA Stand



3 PCS

6.2112.130

Adapter plug contact 5.5 mm / screw contact, gold-plated

For VA electrodes

