



NIRS XDS Process Analyzer – MicroBundle SinglePoint

A629280110C

The NIRS XDS Process Analyzer – MicroBundle provides the next generation of process analyzers for real-time analysis in the pharmaceutical and chemical industries. Non-destructive, accurate measurements of practically any sample type (powders, slurries, liquids,..) are performed directly in the process line, granulator, dryer or reaction vessel. The micro bundle fiber optics allows the measuring point to be up to 35 meter away from the XDS NIR process analyzer.

Phụ kiện tùy chọn

Order no.	Description
A649281190	NIRS XDS Process Contact Lamp 50W, 12 Volt- White lead
	Replacement lamp with contact for NIRS Process instruments. 50W, 12 Volt- White lead.

