



NIRS XDS Process Analyzer – MicroBundle 9 Channels

A629280130C

The NIRS XDS Process Analyzer – MicroBundle provides the next generation of process analyzers for real-time analysis in the pharmaceutical and chemical industries. Non-destructive, accurate measurements of practically any sample type (powders, slurries, liquids,..) are performed directly in the process line, granulator, dryer or reaction vessel. Up to 9 probes and/or flow-cells can be connected to the analyzer. All 9 channels can be configured independently from each other and the micro bundle fiber optics allows the 9 measuring points to be up to 35 meter away from the XDS NIR process analyzer.

Phụ kiện tùy chọn

Order no.	Description
-----------	-------------

A649281190

NIRS XDS Process Contact Lamp 50W, 12 Volt- White lead

Replacement lamp with contact for NIRS Process instruments. 50W, 12 Volt- White lead.

