



NIRS XDS Process Analyzer – SingleFiber SinglePoint

A629280210C

The NIRS XDS Process Analyzer – SingleFiber provides the next generation of process analyzers for real-time analysis in the pharmaceutical and chemical industries. Non-destructive, accurate measurements of clear liquids are performed directly in the process line or in a reaction vessel. Single fiber optics allows the measuring point to be up to 100 meter away from the XDS NIR process analyzer.

Phụ kiện tùy chọn

Order no.	Description
A649281190	NIRS XDS Process Contact Lamp 50W, 12 Volt- White lead
	Replacement lamp with contact for NIRS Process instruments. 50W, 12 Volt- White lead.

