





# NIRS Analyzer PRO – ContactReflection

A629281110

The NIRS Analyzer Pro is a single point process analysis system based on high resolution diode array technology. It provides non-destructive analysis of solids, granules, or powders through a glass window directly in the process line without bypass – a true inline system.

## Phụ kiện tùy chọn

| Order no.  | Description                                                                          |
|------------|--------------------------------------------------------------------------------------|
| 6.7440.480 | NIRS Process, I/O Communication box                                                  |
|            | I/O Communication box for process instruments with 24 analog and 24 digital channels |



A660044277

NIRS Analyzer PRO Lamp 12V/20W

Replacement lamp for the NIRS Analyzer Pro (12V/20W)

