





BWID Standard Raman ID Software

BWT-840000126


Data acquisition, material identification and reporting software for B&W Tek portable Raman spectrometers to rapidly identify or verify materials based on user-built libraries, proprietary B&W Tek and third-party libraries. Includes system performance tests. One user license provided. Free upgrades included within one year of purchase date.

Phụ kiện tùy chọn

Order no.	Description
BWT-840001123	Bluetooth 1D barcode scanner
	Bluetooth 1D barcode scanner for use with BWID software operating a portable Raman instrument



BWID library of illicit materials includes drugs, TICs, TIMs, precursors, and general chemicals.



BWID Metrohm Comprehensive RamanLibrary - More than 27,700 spectral entries of chemicals, minerals, organic chemicals, solvents, polymers, dyes and pigments, TICs, TIMs, inorganics and organometallics, personal care products, and illicit materials.

PTL-MCRL

