





BWID Pharma Raman ID Software

BWT-840000127

Data acquisition and material identification software for B&W Tek portable Raman spectrometers to rapidly identify or verify materials based on user-built libraries, proprietary B&W Tek and third-party libraries. Meets all the requirements of FDA 21 CFR Part 11, including system access security, audit trail of data activities, and system performance tests. One user license provided. Free upgrades included within one year of purchase date.

Phụ kiện tùy chọn

Order no.	Description
BWT-840001123	Bluetooth 1D barcode scanner
	Bluetooth 1D barcode scanner for use with BWID software operating a portable Raman instrument



Metrohm Comprehensive USP Raman Library for BWID containing ~ 900 traceable certified items.

