




96X Screen-Printed Carbon Electrode

DRP-96X110-U4

96X Screen-Printed Carbon Electrode. Formed by 96 three-electrode electrochemical cells. Gold plated contact paths are printed in the backside of the plate

Phụ kiện tùy chọn

Order no.	Description
DRP-CONNECTOR96X	Connects 96X format Screen-Printed Electrodes with any kind of (multi) potentiostat.
	Connects 96X format Screen-Printed Electrodes with any kind of (multi) potentiostat. Two manual rotating selectors allow to select the row letter (from A to H) and the column number (from 1 to 12) identifiers of the wells to be analysed.

