


Cell for Screen-Printed Electrodes in PEEK

DRP-CELL-PEEK

Cell in PEEK for batch analysis which allows optional stirring and it's recommended for deaerated samples and for standard additions mode. Suitable for working with all catalogue Screen-Printed Electrodes but recommended to be used with work in solution format SPEs.

Phụ kiện tùy chọn

Order no.	Description	
DRP-CAC	Cable connector for Screen-Printed Electrodes Connects DropSens Screen-Printed Electrodes to any kind of potentiostat. Recommended when working in batch or flow analysis.	
DRP-CAST	μStat Cable connector for Screen-Printed Electrodes Connects any DropSens equipment with DropSens Screen-Printed Electrodes	
DRP-CELLHOLDER	Cell Holder Accessory to be used in combination with CELL and CELL-PEEK when working in batch analysis for set up improvement.	