



Flow Cell for magnetic assays with Screen-Printed Electrodes

DRP-CFLWCL-MAGN

Methacrylate wall-jet Flow-Cell associated with a powered magnet designed to do FIA with magnetic beads. Suitable to be used with standard format Screen-Printed Electrodes with the electrochemical cell in the middle of the strip. Fittings included.

Phụ kiện tùy chọn

Order no.	Description
-----------	-------------

DRP-110STR-U50	Streptavidin modified Screen-Printed Carbon Electrode
----------------	---


Streptavidin modified Screen-Printed Carbon Electrodes provide a stable high affinity surface for a large amount of biotinylated molecules.



DRP-CAC

Cable connector for Screen-Printed Electrodes

Connects DropSens Screen-Printed Electrodes to any kind of potentiostat. Recommended when working in batch or flow analysis.



DRP-C013-U75

Screen-Printed Silver Electrode (Aux.:C; Ref.:Ag)

Screen-Printed Silver Electrode (Aux.:C; Ref.:Ag). Working Electrode of 1.6 mm diameter.



DRP-CAST

μ Stat Cable connector for Screen-Printed Electrodes

Connects any DropSens equipment with DropSens Screen-Printed Electrodes



DRP-F10-U75

Screen-Printed Ferrocyanide/Carbon Electrode

Screen-Printed Ferrocyanide/Carbon Electrode. Ideal for the determination of hydrogen peroxide at a low detection potential.

