




Flow Cell in PEEK for Screen-Printed Electrodes

DRP-FLWCL-PEEK

PEEK wall-jet Flow-Cell for FIA. Suitable to be used with standard format Screen-Printed Electrodes with the electrochemical cell in the middle of the strip. Closing system with powerful magnets. Fittings included.

Phụ kiện tùy chọn

Order no.	Description
DRP-250BT-U75	Screen-Printed Gold Electrode (Aux.:Pt; Ref.:Ag) / Ink BT
	Screen-Printed Gold Electrode (Aux.:Pt; Ref.:Ag) / Ink BT. Suitable for working with microvolumes, for decentralized assays or to develop specific sensors.



DRP-CAC

Cable connector for Screen-Printed Electrodes

Connects DropSens Screen-Printed Electrodes to any kind of potentiostat. Recommended when working in batch or flow analysis.



DRP-710-U75

Screen-Printed Prussian Blue/Carbon Electrode

Screen-Printed Prussian Blue/Carbon Electrode. Ideal for the determination of hydrogen peroxide at a low detection potential.



DRP-CAST

μStat Cable connector for Screen-Printed Electrodes

Connects any DropSens equipment with DropSens Screen-Printed Electrodes



DRP-C223AT-U75 Screen-Printed Gold Electrode (Aux.:Au; Ref.:Ag) d. 1.6 / Ink AT

Screen-Printed Gold Electrode (Aux.:Au; Ref.:Ag) / Ink AT. Working Electrode of 1.6 mm diameter.

