



797 VA Computrace with automatic standard addition for trace analysis

MVA-02

MVA-02 is a simple to operate, partially automated analysis system for voltammetric trace analysis and for education. For example, heavy metals can be conveniently determined in single samples without a sample changer. The system is comprised of the 797 VA Computrace with two 800 Dosinos for the automatic addition of standard solution and supporting electrolyte.

The system is supplied with extensive accessories for a complete measuring site setup, including VA Computrace software and electrodes: Multi-Mode Electrode pro (MME pro), Ag/AgCl reference electrode and Pt auxiliary electrode.

Phạm vi giao hàng MVA-02

Qt.	Order no.	Description
-----	-----------	-------------

2 PCS

2.800.0010

800 Dosino

The 800 Dosino is a drive with write/read hardware for intelligent Dosing Units. With fixed cable (length 150 cm).



1 PCS

6.1608.050

Drying bottle / 100 mL / GL 45



1 PCS

6.1608.090

Bottle / 250 mL / GL 45



2 PCS

6.1618.060

Thread adapter GL 45 to GL 45

Thread adapter GL45 to GL45; e.g. for reagent bottles with extended flange



1 PCS

6.1805.020

FEP tubing / M6 / 52 cm

With light and kink protection



1 PCS

6.1805.120

FEP tubing / M6 / 100 cm

With light and kink protection



1 PCS

6.1819.010

PTFE tube

As gas inlet in washing bottle of VA Stand



1 PCS

6.2055.100

Bottle holder, double

Bottle holder for 2 x 1 L bottles. For 6.5617.000 aspiration equipment




807 Dosing Unit with integrated data chip with 2 mL glass cylinder and light protection, mountable to a reagent bottle with ISO/DIN GL 45 glass thread. FEP tubing connection, antidiffusion tip.



807 Dosing Unit with integrated data chip with 5 mL glass cylinder and light protection, mountable to a reagent bottle with ISO/DIN GL 45 glass thread. FEP tubing connection, antidiffusion tip.



Phụ kiện tùy chọn

Order no.	Description	
2.843.0040	843 VA Membrane Pump Station for Compact VA systems The 843 Membrane Pump Station has two built-in membrane pumps that can be controlled directly by the 797 VA Computrace via the remote signal. It is used for the automatic rinsing and emptying of the measuring vessel in partly automated and fully automated VA and CVS systems with the 797 VA Computrace. The version for voltammetry is equipped with the complete accessories for the automatic rinsing and emptying of the measuring vessel.	 The image shows the 843 VA Membrane Pump Station, a compact white unit with a green front panel. It is connected via blue tubing to two large white plastic bottles. The bottle on the left is labeled 'Waste' and has a red cap, while the bottle on the right is labeled 'Rinsing Solution' and has a green cap.
2.846.0010	846 Dosing Interface USB-capable control unit for connecting a maximum of four 800 Dosinos or 805 Dosimats for dosing and Liquid Handling tasks. A Touch Control or the connection to a PC with OMNIS Software, tiamo™, MagIC Net, viva, or 797 VA Computrace is required for operation.	 The image shows the 846 Dosing Interface, a white rectangular control unit. It features a green label on the front with the text '846 Dosing Interface' and the Mettler logo. A 'Power' button is visible on the right side of the unit.