



# "ProfIC 12 – Anion" – Professional IC gradient system

ProfIC-12-Anion

The Professional IC system with conductivity detection enables the determination of anions with a binary high-pressure gradient. In addition to step gradients, any linear, convex or concave gradients can be run with this system in order to work in a time-optimized fashion.

Schematic diagram

Phạm vi giao hàng ProfIC-12-Anion

| Qt. | Order no. | Description |
|-----|-----------|-------------|
|-----|-----------|-------------|

Compact and intelligent high performance conductivity detector for intelligent IC instruments. Outstanding temperature stability, the complete signal processing within the protected detector block and the latest generation of DSP – Digital Signal Processing – guarantee the highest precision of the measurement. No change of measuring ranges (not even automatic ones) is required, due to the dynamic working range.



The 858 Professional Sample Processor – Pump processes samples from 500  $\mu\text{L}$  to 500 mL. The sample transfer takes place either with the installed bidirectional two-channel peristaltic pump or with an 800 Dosino.



## Phụ kiện tùy chọn

| Order no.  | Description                                                                                                                                                                                                                                                                                                                                                                |                                                                                       |
|------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
| 6.1014.000 | <p>Metrosep A Trap 1 - 100/4.0</p> <p>High-capacity anion column, which is used to purify the eluent flow. Even chemicals of the highest quality level, e.g. "Ultrapur" or "puriss", can still contain minimal anion contaminations. These are reliably held back by the Metrosep A Trap 1 - 100/4.0.</p> <p>This column is primarily used with gradient applications.</p> |    |
| 6.1608.080 | <p>Bottle / 300 mL</p> <p>Rinsing vessel for IC Sample Processors.</p>                                                                                                                                                                                                                                                                                                     |    |
| 6.2041.440 | <p>Sample rack, 148 x 11 mL and 3 x 300 mL (PVC)</p> <p>Sample rack for 858 Professional Sample Processor, 815 Robotic USB Sample Processor XL, suitable for 148 sample vessels and 3 rinsing bottles. The following sample vessels and rinsing bottles can be used: 6.2743.050, 6.1608.080, 6.1627.000, 6.2743.040.</p> <p>Plastic: Polyvinyl chloride (PVC)</p>          |  |
| 6.2041.760 | <p>Sample rack 54 x 11 mL + 1 x 300 mL</p> <p>Sample rack with small diameter for 54 x 11 mL sample vessels in two rows and additionally 1 x 300 mL bottle.</p>                                                                                                                                                                                                            |  |

6.2743.040 Vial PP, 2.5 mL, 2000 pieces

2.5 mL. 2000 pieces. For IC Sample Processors.

The IC sample vessels are made of polypropylene (PP) plastic for health-care applications. Production and packaging both take place under sterile and clean-room conditions (Clean Room Class 7, in accordance with DIN ISO 14644). This guarantees minimal degrees of contamination. The sample vessels are subjected to random checks for leachable anions and cations.



---

6.2743.050 Sample vessel 11 mL

2000 pieces. For IC Sample Processors and VA Autosamplers.

The IC sample vessels are made of polypropylene (PP) plastic for health-care applications. Production and packaging both take place under sterile and clean room conditions (Clean Room Class 7, in accordance with DIN ISO 14644). This guarantees minimal degrees of contamination. The sample vessels are subjected to random checks for leachable anions and cations.



---

6.2743.070 Stopper with perforation

2000 pieces. For IC Sample Processors



---

6.5333.000 IQ/OQ Kit for IC

The IQ/OQ Kit contains all parts and standard solutions required for IQ /OQ in ion chromatography.

