



NIRS XDS Interactance OptiProbe 模块

2.921.0510

坚固耐用的测量系统可用于监测工业设备和实验场所中的反应

NIRS XDS Interactance OptiProbe Analyzer 分析仪可方便可靠地监控工业设备和实验场所中的化学反应。在方法开发以及生产过程的工艺放大等应用领域中，NIRS XDS Interactance OptiProbe 分析仪可提供不同样品的识别和品质测定的精确结果。在使用反射探针测量固体、高度散射液体和浆料的过程中，分析水性制品、透明液体和溶剂时可使用透射反射探针。探针通过光纤与分析仪相连，这样即使在恶劣或危险的过程环境条件下也可确保可靠地测量。

Scope of delivery 2.921.0510

Qt.	Order no.	Description
-----	-----------	-------------

坚固耐用的测量系统可用于监测工业设备和实验场所中的反应

NIRS XDS Interactance OptiProbe Analyzer 分析仪可方便可靠地监控工业设备和实验场所中的化学反应。在方法开发以及生产过程的工艺放大等应用领域中，NIRS XDS Interactance OptiProbe 分析仪可提供不同样品的识别和品质测定的精确结果。在使用反射探针测量固体、高度散射液体和浆料的过程中，分析水性制品、透明液体和溶剂时可使用透射反射探针。探针通过光纤与分析仪相连，这样即使在恶劣或危险的过程环境条件下也可确保可靠地测量。



Optional accessories

Order no.	Description
6.7430.080	XDS 可更换样品套件 可更换探针组件包括一个反射探针、一个透射探针和各种间距支架，配有固定螺栓用于调节不同路径长度（2 cm、4 cm、5 cm 和 10 cm）。此外还有一个探针保护罩、一把内六角扳手和一个密封环。 这套组件是 XDS SmartProbe Analyzer 分析仪和 XDS Interactance Optiprobe Analyzer 分析仪的组成部分。 The image shows the components of the XDS replaceable probe kit. It includes two long, thin, cylindrical metal probes, one of which is shown with a protective cap. There are also four small, cylindrical metal spacers of different lengths, and a small hex key (Allen key) and a circular seal ring.