



823 Membrane Pump Unit

2.823.0010

Membranpumpe zum direkten Anschluss an Sample Processors.

Lieferumfang 2.823.0010

Qt.	Order no.	Beschreibung
1 STK	1.823.0010	823 Membrane Pump Unit Membranpumpe zum direkten Anschluss an Sample Processors.



2 STK

6.1801.120

PVC-Schlauch / 4 mm / 6 mm / 2 m



1 STK

6.1812.000

PTFE-Schlauch 4 mm / 6 mm / 4 m



1 STK

6.1828.000

PVDF-Anschlussnippel

Für Kanister 6.1621.000.

