



Eco Dosimat

2.1007.0010

The Eco Dosimat is a dispensing unit with built-in magnetic stirrer and touch display for universal use.

Increase your efficiency and precision for classic laboratory tasks such as

- the precise dosing of any volume.
- manual titrations.
- the preparation of ion standards.
- controlled dosing with specified dosing rate.

A report in accordance with GLP can be stored as a PDF or printed out for all of the activities listed.

Including push-button cable for manual titration and dispensing control.

Freely selectable cylinder unit with 5, 10, 20, and 50 mL.

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 2.1007.0010

Qt.	Order no.	Description
-----	-----------	-------------

1 PCS	1.1007.0010	Eco Dosimat
-------	-------------	-------------

The Eco Dosimat is a dispensing unit with built-in magnetic stirrer and touch display for universal use.

Increase your efficiency and precision for classic laboratory tasks such as

- the precise measurement of volumes.
- of manual titration.
- the creation of ion standards.
- controlled dosing with specified dosing rate.

A report in accordance with GLP can be stored as a PDF or printed out for all of the activities listed.

Including push-button cable for manual titration and dispensing control.

Freely selectable cylinder unit with 5, 10, 20, and 50 mL.



1 PCS	6.1446.030	Ball stopper / B-14/(15)
-------	------------	--------------------------

For the movable holding of buret tips in titration vessels



1 PCS	6.1446.080	Threaded stopper / M8
-------	------------	-----------------------



1 PCS

6.1543.060

Tip / M6 thread

Tip with M6 thread. Together with the 6.1805.160 capillary and the 6.1446.030 link stopper, it makes up the entire 6.1537.010 buret tip.



1 PCS

6.1543.200

Titration tip / M6 thread

Titration tip with antidiffusion valve.



1 PCS

6.1602.105

Bottle attachment / GL 45 / green

For direct dosing from reagent bottles with GL 45 thread.



1 PCS

6.1619.010

Adsorber complete, for Exchange Unit



2 PCS

6.1805.090

FEP tubing / M6 / 31 cm

With light and kink protection



1 PCS	6.1805.100	FEP tubing / M6 / 40 cm
		With light and kink protection



1 PCS	6.1819.020	Tube
		For M6 threaded holes. For 6.1602.1XX bottle attachments



1 PCS	6.1903.020	Stirring bar / 16 mm
		Stirring Bar with magnetic core, PTFE coated, length 16 mm



1 PCS	6.1903.030	Stirring bar / 25 mm
		Stirring bar with magnetic core, PTFE covering, length 25 mm.



1 PCS

6.2013.010

Clamping ring

For support rods with a diameter of 10 mm.



1 PCS

6.2016.070

Support rod / 400 mm



1 PCS

6.2021.020

Electrode holder

Electrode holder for 4 electrodes and 2 buret tips



1 PCS

6.2023.020

Clip for SGJ 14/15

Clip for SGJ 14/16



1 PCS

6.2709.070

Guiding sleeve

For fixation of buret and dosing tips in titration heads and electrode holders



1 PCS	6.2107.100	Manual Dosing Controller
-------	------------	--------------------------

Handle for manual dosing and manual titration with the Eco Dosimat or Dosimat plus



1 PCS	6.1546.040	Piston tool for Eco Dosimat/Eco Titrator
-------	------------	--

Tool for removing and mounting the piston.



1 PCS	6.2164.010	Power supply unit 100 - 240 V/24 V DC
-------	------------	---------------------------------------

Power supply unit for 915 KF Ti-Touch, 916 Ti-Touch, 917 Coulometer and all Eco Titrators / Eco Dosimat.



1 PCS	6.1608.023	Amber glass bottle / 1000 mL / GL 45
-------	------------	--------------------------------------

For Exchange Units. Bottle for auxiliary solutions.



Optional accessories

Order no.	Description
6.03002.150	Cylinder unit Eco 5 mL Cylinder unit 5 mL for Eco Dosimat / Eco Titrator. Without dosing tubing and antidiffusion tip.
6.03002.210	Cylinder unit Eco 10 mL Cylinder unit 10 mL for Eco Dosimat / Eco Titrator. Without dosing tubing and antidiffusion tip.
6.03002.250	Cylinder unit Eco 50 mL Cylinder unit 50 mL for Eco Dosimat / Eco Titrator. Without dosing tubing and antidiffusion tip.
6.03002.220	Cylinder unit Eco 20 mL Cylinder unit 20 mL for Eco Dosimat / Eco Titrator. Without dosing tubing and antidiffusion tip.
2.142.0100	Custom Q3X thermal printer



Compact printer with USB interface for

- 900 Touch Control
- 915 KF Ti-Touch
- 916 Ti-Touch
- 917 Coulometer
- 877 / 848 Titrino plus
- 865 / 876 Dosimat plus
- 91X Meter (cable 6.2151.140)
- Eco Dosimat / Titrator
- 862 Compact Titrosampler
- 870 KF Titrino plus
- 899 Coulometer



Paper width 60 mm (40 characters). Including 6.2151.120 USB cable.

2.143.0200

TX-900MH impact printer

Compact impact printer with USB and RS232 interface for

- 900 Touch Control (+)
- 915 KF Ti-Touch (+)
- 916 Ti-Touch (+)
- 917 Coulometer (+)
- 91X Meter (+)
Eco Dosimat/Titrator (+)
- 877/848 Titrino plus (additional 6.2151.100 required)
- 865/876 Dosimat plus (additional 6.2151.100 required)
- 862 Compact Titrosampler (additional 6.2151.100 required)
- 870 KF Titrino plus (additional 6.2151.100 required)
- 899 Coulometer (additional 6.2151.100 required)



(+) Cable included in scope of delivery
