



OMNIS Oven Module 2R

2.1017.0010

Module for installation in the module plate of the OMNIS Sample Robot Oven. This workstation picks up the 2R sample vial and transports it from the handover position into the oven and then back again after the analysis is completed.

Below, the accessories are grouped into Scope of delivery and Optional accessories. Please keep this printout at hand for ordering replacement material. These lists may be subject to change.

Scope of delivery 2.1017.0010

Qt.	Order no.	Description
1 PCS	1.1017.0010	OMNIS Oven Module 2R
		OMNIS Oven Module 2R



Optional accessories

Order no.	Description
6.02041.140	OMNIS sample rack 77 x 2R OMNIS sample rack for OMNIS Sample Robot Oven, suitable for 77 2R sample vials.

