

OMNIS Sample Robot Oven

2.1030.0010

Main module for setting up an OMNIS Sample Robot Oven in size S. This module is equipped with the module and rack plate. It is already equipped with the main lift and a gripper. It must be equipped with sample racks and at least 1 oven module in order to use it to set up a functioning Sample Robot Oven.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Optional accessories

Order no.	Description
6.02041.140	OMNIS sample rack 77 x 2R
	OMNIS sample rack for OMNIS Sample Robot Oven, suitable for 77 2R sample vials.

