



Ionizer for 864 Robotic Balance Sample Processor, 230V

2.136.0306

Ionizer (HAUG) for 864 Robotic Balance Sample Processor with power supply for 230V.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Scope of delivery 2.136.0306

Qt.	Order no.	Description
-----	-----------	-------------

1 PCS

1.136.0306

Ionizer for 864 Robotic Balance Sample Processor,
230V

Ionizer (HAUG) for 864 Robotic Balance Sample Processor with power
supply for 230V.



1 PCS

6.9013.000

Ionisation rod for ionizer

The ionisation rod is mounted with the holder (6.2057.160) to the 864
Robotic Balance Sample Processor.



Optional accessories

Order no.	Description	
6.2057.160	Holder for ionization rod for 864 Robotic Balance Sample Processor	
