



USB thermal printer Neo's

2.141.0100

Compact printer with USB interface for the Ti-Touch, the Titrino plus and Dosimat plus and the 91X Meter instruments.

Paper width 60 mm (40 characters). Including 6.2151.120 USB cable.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Scope of delivery 2.141.0100

Qt.	Order no.	Description
-----	-----------	-------------

1 PCS

1.141.0100

USB thermal printer Neo's

Compact printer with USB interface for the Ti-Touch, the Titrino plus and Dosimat plus and the 91X Meter instruments.

Paper width 60 mm (40 characters). Including 6.2151.120 USB cable.



1 PCS

6.2151.120

Cable USB MINI (OTG) - USB B, 1.8 m

For connecting USB instruments.



Optional accessories

Order no.	Description	
6.9001.050	Thermal paper for Custom Neos-N / Q3X	
	Thermal paper 60 mm / 25 m for Custom Neos-N / Q3X. New: 25-year legibility with correct storage.	
	1 roll	

