



Custom Q3X thermal printer

2.142.0100

Compact printer with USB interface for

- 900 Touch Control
- 915 KF Ti-Touch
- 916 Ti-Touch
- 917 Coulometer
- 877 / 848 Titrino plus
- 865 / 876 Dosimat plus
- 91X Meter (cable 6.2151.140)
- Eco Dosimat / Titrator
- 862 Compact Titrosampler
- 870 KF Titrino plus
- 899 Coulometer

Paper width 60 mm (40 characters). Including 6.2151.120 USB cable.

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 2.142.0100

Qt.	Order no.	Description
1 PCS	1.142.0100	Custom Q3X thermal printer
		Compact printer with USB interface for
		• 900 Touch Control • 915 KF Ti-Touch • 916 Ti-Touch • 917 Coulometer • 877 / 848 Titrino plus • 865 / 876 Dosimat plus • 91X Meter • Eco Dosimat / Titrator
		Paper width 60 mm (40 characters). Including 6.2151.120 USB cable.
1 PCS	6.2151.120	Cable USB MINI (OTG) - USB B, 1.8 m
		For connecting USB instruments.
1 PCS	6.9001.050	Thermal paper for Custom Neos-N / Q3X
		Thermal paper 60 mm / 25 m for Custom Neos-N / Q3X. New: 25-year legibility with correct storage.
		1 roll



Optional accessories

Order no.	Description	
6.9001.050	Thermal paper for Custom Neos-N / Q3X Thermal paper 60 mm / 25 m for Custom Neos-N / Q3X. New: 25-year legibility with correct storage. 1 roll	
6.2151.140	Y cable USB A St - USB B St - Mini B St Y cable for connection of a USB printer to the pH/Conductometers 912 / 913 / 914. This cable allows a printer and the power supply unit to be connected to the measuring instrument at the same time.	