




Sample changer 760

2.760.0010

The sample changer for simple automation tasks that can be solved with standard programs. Sample vessels and sample racks are identical to those with the Sample changer 730. Pumps can be actuated externally if required.

Below, the accessories are grouped into Scope of delivery and Optional accessories. Please keep this printout at hand for ordering replacement material. These lists may be subject to change.

Scope of delivery 2.760.0010

Qt.	Order no.	Description
1 PCS	6.2752.010	Plug cover
		Used with the 717, 730, 760 Sample Changers, IC Sample Processors, 774 Oven Sample Processor



2 PCS

6.1236.020

Sleeve with SGJ 14/12 mm

Sleeve with SGJ 14/12 mm and O-ring.



1 PCS

6.1446.160

Thermometer stopper, B-14/15, oblique insertion

To have the inserts of the 6.1458.010 titration head in a sloping position.



1 PCS

6.1805.110

FEP tubing / M6 / 80 cm

With light and kink protection



1 PCS

6.1815.010

Spiral band 0.5 m

For holding together different cables or tubing.



1 PCS

6.2041.310

Sample rack, 12 x 250 mL, (PVC)

Sample rack for 810 Sample Processor, 814 USB Sample Processor, 855 Robotic Titrosampler "Light", suitable for 12 sample beakers. The following sample beakers can be used: 6.1432.320, 6.1432.323, 6.1453.220, 6.1453.250.

Plastic: Polyvinyl chloride (PVC)



3 PCS

6.2709.070

Guiding sleeve

For fixation of buret and dosing tips in titration heads and electrode holders




1 PCS

6.5610.010

KFT equipment for 810/814 Sample Processor

Karl Fischer equipment for 810 or 814 Sample Processor.



Optional accessories

Order no.	Description	
6.1453.220	Sample beaker, plastic (PP), 200 mL, 1 piece Sample beaker, plastic (polypropylene, PP), 200 mL, suitable for sample racks: 6.2041.310, 6.2041.400, 6.2041.430, 6.2041.510, 6.2041.820, 6.2068.020, and 6.2068.030, 1 piece	
6.1453.250	Sample beaker, plastic (PP), 250 mL, 1 piece Sample beaker, plastic (polypropylene, PP), 250 mL, suitable for sample racks: 6.2041.310, 6.2041.400, 6.2041.430, 6.2041.510, 6.2041.820, 6.2068.020, and 6.2068.030, 1 piece	
6.1458.010	Titration head, 6 x SGJ 14 and 3 x SGJ 9 openings. With 6xSGJ 14 and 3x SGJ 9 openings.	
6.1458.020	Titration head / 4 x M10 connectors With 4 M10 screw nipples, including stoppers.	

6.1458.030

KF titration head /4 x M10 connectors

With 4 M10 screw nipples, including stoppers.