




Sample changer KFT 760

2.760.0020

Complete system for automatic volumetric Karl Fischer titration.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Scope of delivery 2.760.0020

Qt.	Order no.	Description
1 PCS	6.2752.010	Plug cover
		Used with the 717, 730, 760 Sample Changers, IC Sample Processors, 774 Oven Sample Processor




1 PCS

1.741.0010

741 Magnetic Stirrer

Magnetic stirrer for sample changers.

The magnetic stirrer is located beneath the sample rack.



1 PCS

6.0340.000

Double Pt-wire electrode for sample changer (fixed cable)

Indicator electrode (diameter = 5.3 mm, fixed cable 2 m) used for automated Karl Fischer titration.



24 PCS

6.1432.210

Sample beaker, clear glass, 75 mL, 1 piece

Sample beaker, clear glass, 75 mL, suitable for sample racks: 6.2041.340, 6.2041.350, 6.2041.800 and 6.2041.900, 1 piece



1 PCS

6.1458.030

KF titration head / 4 x M10 connectors

With 4 M10 screw nipples, including stoppers.



2 PCS 6.1805.110 FEP tubing / M6 / 80 cm
With light and kink protection



1 PCS 6.1815.010 Spiral band 0.5 m
For holding together different cables or tubing.



1 PCS 6.1821.000 Aspiration tube
For KF Sampler



24 PCS 6.1903.010 Stirring bar / 12 mm
Stirring Bar with magnetic core, PTFE coated, length 12 mm



24 PCS

6.2037.040

Foil holder

Foil holder for 75 mL sample beaker for Karl Fischer for KF Sample Changer



1 PCS

6.2041.340

Sample rack, 24 x 75 mL, (PVC)

Sample rack for 810 Sample Processor, 814 USB Sample Processor, 855 Robotic Titrosampler "Light", suitable for 24 sample beakers. The following sample beakers can be used: 6.1432.210, 6.1459.400.

Plastic: Polyvinyl chloride (PVC)



1 PCS

6.2820.000

Al foil

1000 circular sheets of 80 mm diameter. Thickness: 0.010 mm. Used together with 6.2037.030 Foil holder for 6.1432.210 Sample beaker. For 674 KF Sample Changer

