

# Pump Unit 772 "Aspiration"

2.772.0020

Peristaltic pump for the aspiration of titrated solution and suspensions. The "Aspiration" version contains extensive accessories for aspiration. The pump can be actuated through a Relay Box 731. It is suitable for the conveyance of solvents.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

## Scope of delivery 2.772.0020


| Qt.   | Order no.  | Description                                                                                                                    |                                                                                       |
|-------|------------|--------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
| 1 PCS | 6.1446.160 | Thermometer stopper, B-14/15, oblique insertion<br>To have the inserts of the 6.1458.010 titration head in a sloping position. |   |
| 1 PCS | 6.1543.170 | Aspiration tip / M8 thread<br>Aspiration tip for the Sample Processors.                                                        |  |
| 1 PCS | 6.1621.000 | Container 10 L<br>As rinsing or waste container in automated systems.                                                          |  |

1 PCS

6.1805.510

PTFE tubing / M8 / 60 cm

With protection against kinks.



---

1 PCS

6.1812.000

PTFE tubing 4 mm / 6 mm / 4 m



---

2 PCS

6.1820.050

M6 / Olive adapter

For 6.1826.100 pump tubing, to peristaltic pumps.



---

1 PCS

6.1826.100

Pump tubing 6.4 mm / 1 m

Replaces the pump tubing used with the 772 Pump Unit and the 843 Pump Station (peristaltic).



1 PCS

6.1828.000

PVDF connection nipple

For 6.1621.000 canister.



---

1 PCS

6.2621.140

Hexagon key 2.5 mm



# Optional accessories

| Order no.  | Description            |                                                                                    |
|------------|------------------------|------------------------------------------------------------------------------------|
| 6.1826.160 | Viton pump tubing      |                                                                                    |
|            | For peristaltic pumps. |  |