



786 Swing Head with titration swing arm, left- or right-swinging

2.786.0030

Additional equipment for Robotic Sample Processors XL for direct titration in sample vessels of 75 mL and more.

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 2.786.0030

Qt.	Order no.	Description
-----	-----------	-------------

1 PCS

1.786.0010

786 Swing Head with transfer head left

Additional equipment for Robotic Sample Processors XL for transferring / pipetting from small sample vessels into larger titration vessels on the rack or into an external titration cell; this can be located to the left of the work station.



1 PCS

6.1462.050

Robotic arm with titration head for 786 Swing Head,
left or right swinging

Titration head for Swing Head



1 PCS

6.1808.040

Thread adapter M6 outer / M8 inner

Outer thread M6, inner thread M8.



1 PCS

6.1909.030

Stirring propeller / 95 mm

Propeller stirrer with PP shaft, PVDF tip. Length from the SGJ: 104 mm. For KF applications with 722, 802 Rod Stirrer



2 PCS

6.2042.030

Buret tip clip (rod stirrer)

Clip for proper alignment of buret tips, dosing or aspiration tubings using a rod stirrer in the titration beaker of a sample changer



1 PCS

6.2621.120

Hexagon key 1.5 mm

