



# 804 Ti Stand without stand rod

2.804.0010

Titration stand and controller for 802 Stirrer. The 804 Ti Stand together with the optional 802 Stirrer provides an alternative to the magnetic stirrer. Without stand.

Below, the accessories are grouped into Scope of delivery and Optional accessories.  
Please keep this printout at hand for ordering replacement material.  
These lists may be subject to change.

## Scope of delivery 2.804.0010

| Qt.   | Order no.  | Description                                                                                                                                                             |
|-------|------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1 PCS | 1.804.0010 | 804 Ti Stand without stand rod                                                                                                                                          |
|       |            | Titration stand and controller for 802 Stirrer. The 804 Ti Stand together with the optional 802 Stirrer provides an alternative to the magnetic stirrer. Without stand. |



# Optional accessories

| Order no.  | Description                                                                                                                          |                                                                                     |
|------------|--------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
| 6.2001.090 | Base plate on both sides<br><br>Double stand for attaching 2 x 801, 803, 804 to a Titrando                                           |  |
| 6.2151.010 | Cable MSB plug / socket<br><br>Extension cable for MSB connections and for Touch Control - Titrando, Mini-DIN socket - mini-DIN plug |  |