



823 Membrane Pump Unit

2.823.0010

Membrane pump for direct connection to Sample Processors.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Scope of delivery 2.823.0010

Qt.	Order no.	Description
1 PCS	1.823.0010	823 Membrane Pump Unit
		Membrane pump for direct connection to Sample Processors.



2 PCS

6.1801.120

PVC tubing / 4 mm / 6 mm / 2 m



1 PCS

6.1812.000

PTFE tubing 4 mm / 6 mm / 4 m



1 PCS

6.1828.000

PVDF connection nipple

For 6.1621.000 canister.

