



826 pH mobile with carrying case

2.826.0110

Portable pH meter with IrDa interface. With carrying case, ample accessories for calibration, including 6.0228.020 Primatrode.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Scope of delivery 2.826.0110

Qt.	Order no.	Description
1 PCS	1.826.0010	826 pH mobile
		Portable pH meter with IrDa interface, including carrying strap and batteries.



1 PCS 6.2050.000 Carrying strap for pH Meter
Carrying strap for 704, 826 pH meter.



1 PCS 6.2716.040 Case for 826 pH Meter
Case for 826 pH Meter and accessories



1 PCS 6.1236.050 Sleeve with SGJ 14/12 mm
Sleeve with SGJ 14/12 mm, polyethylene.



1 PCS 6.1446.000 SGJ stopper / B-14/(15)



1 PCS

6.1613.020

Bottle / 25 mL / label pH 7 (empty)



1 PCS

6.1613.030

Bottle / 25 mL / label pH 4 (empty)



1 PCS

6.1614.000

Wash bottle / 250 mL



3 PCS

6.1614.030

Plastic bottle / 50 mL

Plastic bottle



1 PCS

6.2133.000

Batteries LR6, 1.5V, 4 pieces

1.5 V LR6 type. Set of 4 pieces. Used with 704 pH Meter, 773 pH/mV Simulator



1 PCS

6.2307.230

Buffer solutions pH 4, 7 and 9

Mixed buffer solutions pH 4.00/7.00/9.00 (25 °C) in single use sachets, colourless, box of 3 x 10 x 30 mL



1 PCS

6.2308.050

Electrolyte 3 mol/L KCl (50 mL)

Electrolyte solution $c(\text{KCl}) = 3 \text{ mol/L}$ (for Ag/AgCl reference systems)



1 PCS

6.2621.140

Hexagon key 2.5 mm



1 PCS

6.2717.000

PP beaker, 100 mL

Beaker made of PP, 100 mL.



Optional accessories

Order no.	Description
6.2325.000	pHit kit Maintenance kit for pH electrodes The kit contains: • 50 mL cleaning solution • 50 mL 3M KCl solution • 50 mL storage solution • 2 Storage vessels • Instructions for use

