

pH Lab Irda 120 V US

2.827.0015

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 2.827.0015

Qt.	Order no.	Description	
1 PCS	6.2161.020	Plug-in power supply unit US 120V/6V DC Power supply for 827 pH lab	
1 PCS	6.2013.010	Clamping ring For support rods with a diameter of 10 mm.	
1 PCS	6.2016.050	Support rod / 300 mm	

1 PCS

6.2021.020

Electrode holder

Electrode holder for 4 electrodes and 2 buret tips



1 PCS

6.2621.070

5 mm hex key for IC Sample Processors

