



833 Advanced IC Liquid Handling Dialysis Unit

2.833.0040

Dialysis sample preparation module for the determination of anions and cations in difficult matrices such as milk, fruit juice, emulsions, suspensions, etc. Including 2-channel peristaltic pump. Cellulose acetate membrane 0.25 µm.

Below, the accessories are grouped into Scope of delivery and Optional accessories. Please keep this printout at hand for ordering replacement material. These lists may be subject to change.

Scope of delivery 2.833.0040

Qt.	Order no.	Description
7 PCS	6.1803.040	PTFE capillary 0.5 mm i.d. / 1 m Capillary for sample handling in IC.



1 PCS

6.1831.050

PEEK capillary 0.5 mm i.d., 40 cm

For inline dialysis and for sample changers. For IC Dialysis Unit, IC Sample Processor, Compact IC Autosampler, IC Dialysis Sample Processor, IC Liquid Handling Dialysis Unit, Compact VA Autosampler



1 PCS

6.1831.060

PEEK capillary 0.5 mm i.d., 1 m

For inline dialysis and for Sample Changers. For IC Dialysis Unit, IC Sample Processor, Compact IC Autosampler, IC Filtration Sample Processor, IC Dialysis Sample Processor, IC Liquid Handling Dialysis Unit



1 PCS

6.2057.010

Holder for dialysis cell

For the dialysis cell in IC column chambers



1 PCS

6.2128.180

Cable 762/830 - 752/753/754/833

Connecting IC Interface - IC Pump Unit, IC Suppressor Module, IC Dialysis Unit or IC Liquid Handling Unit





Nominal pore size : 0.2 µm. Set of 50 pieces. For inline dialysis



With 4 x 6.2744.060 Blind stoppers. For inline dialysis



With UNF 10/32 connection. For IC instruments, inline dialysis (connection of PTFE capillaries)



2 PCS

6.2744.010

Pressure screw 5x

With UNF 10/32 connection. For the connection of PEEK capillaries



2 PCS

6.2744.034

Coupling olive/UNF 10/32, 2x

Connection between pressure screw and pump tubing. 2 pieces. For IC instruments with peristaltic pump.




4 PCS

6.2744.160

Nozzle/UNF 10/32 coupling with security device

Security device for connecting pump tubing to a nozzle.



1 PCS


6.2831.010

Plastic tweezers

For dialysis and filtration membranes



Optional accessories

Order no.	Description	
6.1826.060	Pump tubing (orange/yellow) For Compact IC , Personal IC , Basic IC	
6.2714.030	Dialysis membrane (polyamide) For inline dialysis and filtration	
6.2755.000	Tubing cartridge For the peristaltic pumps of IC instruments. Contact pressure adjustable.	