



833 IC Liquid Handling Ultra Filtration Unit

2.833.0050

For inline ultrafiltration. Guarantees the reliable removal of particles from the sample. 2-channel peristaltic pump and separate ultrafiltration cell, for example for installation in the 820 IC Separation Center.

Below, the accessories are grouped into Scope of delivery and Optional accessories. Please keep this printout at hand for ordering replacement material. These lists may be subject to change.

Scope of delivery 2.833.0050

Qt.	Order no.	Description
1 PCS	6.1803.040	PTFE capillary 0.5 mm i.d. / 1 m Capillary for sample handling in IC.



1 PCS

6.1803.050

PTFE capillary 0.5 mm i.d. / 20 cm

Capillary for inline dialysis and filtration, for Dialysis Unit, IC Filtration Sample Processor



2 PCS

6.1803.070

PTFE capillary 0.97 mm i.d. / 40 cm

Capillary for filtration, for IC Filtration Sample Processor, IC Liquid Handling Filtration Unit



2 PCS

6.1803.080

PTFE capillary 0.97 mm i.d. / 1 m

Capillary for filtration, for IC Filtration Sample Processor, IC Liquid Handling Filtration Unit



1 PCS

6.1803.100

PTFE capillary, 0.5 mm ID / 40 cm

Capillary for filtration and dialysis



1 PCS

6.2128.180

Cable 762/830 - 752/753/754/833

Connecting IC Interface - IC Pump Unit, IC Suppressor Module, IC Dialysis Unit or IC Liquid Handling Unit



1 PCS

6.2621.070

5 mm hex key for IC Sample Processors



1 PCS

6.2714.020

Filtration membrane reg. cellulose, 0.2 µm, 50 pieces

Nominal pore size: 0.2 µm. Set of 50 pieces. For inline filtration.



1 PCS

6.2744.000

Pressure screw PVDF

With UNF 10/32 connection. For IC instruments, inline dialysis (connection of PTFE capillaries)



1 PCS

6.2744.010

Pressure screw 5x

With UNF 10/32 connection. For the connection of PEEK capillaries



1 PCS

6.2744.030

Coupling olive/UNF 10/32, 4x

Connection between pressure screw and pump tubing. 4 pieces. For IC instruments with peristaltic pump



1 PCS


6.2831.010

Plastic tweezers

For dialysis and filtration membranes



Optional accessories

Order no.	Description	
6.1803.030	PTFE capillary 0.5 mm i.d. / 3 m Capillary for inline dialysis, for Dialysis Unit, IC Dialysis Sample Processor, IC Liquid Handling Dialysis Unit	
6.1826.060	Pump tubing (orange/yellow) For Compact IC , Personal IC , Basic IC	
6.2755.000	Tubing cartridge For the peristaltic pumps of IC instruments. Contact pressure adjustable.	