



843 Pump Station (peristaltic)

2.843.0120

The 843 Pump Station (peristaltic) has two built-in peristaltic pumps. These can be controlled directly via the interface using remote signals or manually by pressing a button.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Scope of delivery 2.843.0120

Qt.	Order no.	Description
1 PCS	1.843.0120	843 Pump Station (peristaltic)
		The 843 Pump Station (peristaltic) has two built-in peristaltic pumps. These can be controlled directly via the interface using remote signals or manually by pressing a button.



4 PCS

6.1801.120

PVC tubing / 4 mm / 6 mm / 2 m



2 PCS

6.1812.000

PTFE tubing 4 mm / 6 mm / 4 m



4 PCS

6.1820.050

M6 / Olive adapter

For 6.1826.100 pump tubing, to peristaltic pumps.



1 PCS

6.1826.100

Pump tubing 6.4 mm / 1 m

Replaces the pump tubing used with the 772 Pump Unit and the 843 Pump Station (peristaltic).



2 PCS

6.1828.000

PVDF connection nipple

For 6.1621.000 canister.



Optional accessories

Order no.	Description	
6.2141.260	Remote cable Compact Sample Changer - Dosimat plus and Titrino plus Cable for remote connection between Compact Sample Changer and Dosimat plus and further devices (e.g. Titrino plus or Pump Station).	
6.2141.300	Remote cable 843 Pump Station to Sample Processor Remote cable for the direct connection of the 843 Pump Station to the pump connectors of the Sample Processors	
6.2761.000	Splash protection for 843 Pump Station (Peristaltic) The peristaltic pump runs with mounted splash protection only.	