



847 USB Lab Link

2.847.0010

847 USB Lab Link - rapid network connection for the Titrande (only for 840 Touch Control).

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Scope of delivery 2.847.0010

Qt.	Order no.	Description
1 PCS	1.847.0010	847 USB Lab Link
		847 USB Lab Link - rapid network connection for the Titrande (only for 840 Touch Control).



1 PCS

6.2151.030

Cable USB A - USB B / 0.6 m

Connecting cable USB A - USB B

