



849 Level Control metal-free (for canisters)

2.849.0020

Additional equipment for Metrohm instruments. For monitoring the filling level of bottles for eluent, auxiliary or waste solutions via remote control. The 849 Level Control prevents pumps from running dry and/or canisters from overflowing and is suitable for use with aqueous solutions, solvents and suspensions.

The metal-free sensor works in solutions with very low conductivity. Therefore it may be used in ion chromatography.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Scope of delivery 2.849.0020

Qt.	Order no.	Description
-----	-----------	-------------

1 PCS

1.849.0020

849 Level Control metal-free (for canisters)

Additional equipment for Metrohm instruments. For monitoring the filling level of bottles for eluent, auxiliary or waste solutions via remote control. The 849 Level Control prevents pumps from running dry and/or canisters from overflowing and is suitable for use with aqueous solutions, solvents and suspensions.

The metal-free sensor works in solutions with very low conductivity. Therefore it may be used in ion chromatography.




1 PCS

6.1113.100

Level sensor MF for 10 L canister

Sensor for Level Control MF



1 PCS

6.1619.000

Adsorber tube for 807 Dosing Unit



1 PCS

6.2151.060

Cable for Level Control and Eluent Production Module

Connecting cable for sensor and Level Control or Eluent Production Module.

