



849 Level Control metal-free (for bottles)

2.849.0030

Additional equipment for automatic "Inline Eluent Preparation" and for remotely monitoring the fill level in eluent or auxiliary solution bottles. The 849 Level Control prevents pump from running dry or prevents it from running over with the optionally available "FULL Sensor". The enclosed "EMPTY sensor" is suitable for use with aqueous solutions, solvents, and suspensions; the length of the rods is intended for the 2 L eluent bottle.

The metal-free sensor works in solutions with very low conductivity. It is therefore suitable for use in ion chromatography, particularly for ultrapure water as well.

Below, the accessories are grouped into Scope of delivery and Optional accessories. Please keep this printout at hand for ordering replacement material. These lists may be subject to change.

Scope of delivery 2.849.0030

Qt.	Order no.	Description
-----	-----------	-------------

1 PCS

1.849.0020

849 Level Control metal-free (for canisters)

Additional equipment for Metrohm instruments. For monitoring the filling level of bottles for eluent, auxiliary or waste solutions via remote control. The 849 Level Control prevents pumps from running dry and/or canisters from overflowing and is suitable for use with aqueous solutions, solvents and suspensions.

The metal-free sensor works in solutions with very low conductivity. Therefore it may be used in ion chromatography.



1 PCS

6.1619.000

Adsorber tube for 807 Dosing Unit



1 PCS

6.2151.060

Cable for Level Control and Eluent Production Module

Connecting cable for sensor and Level Control or Eluent Production Module.



1 PCS

6.1626.000

Bottle attachment GL 45 for level sensor

This bottle attachment, fitted with the 6.2769.XXX sensors, results in the level sensors EMPTY and FULL, respectively, for bottles and canisters.



1 PCS

6.1815.010

Spiral band 0.5 m

For holding together different cables or tubing.



1 PCS

6.2769.110

EMPTY sensor for level sensor for 2 L bottles

This sensor, mounted in the 6.1626.000 bottle attachment, results in the level sensor EMPTY for 2 L bottles. Included are two sensor rods, one aspiration tubing and two holders for fixing the tubing to one of the sensor rods.



Optional accessories

Order no.	Description	
6.2769.000	FULL sensor for level sensor	This sensor, mounted in the 6.1626.000 bottle attachment, results in the level sensor FULL for 1 L and 2 L bottles and for canisters. Included are two sensor rods, one inlet tubing and one holder for fixing the tubing to one of the sensor rods. 
6.2769.100	EMPTY sensor for level sensor for 1 L bottles	This sensor, mounted in the 6.1626.000 bottle attachment, results in the level sensor EMPTY for 1 L bottles. Included are two sensor rods, one aspiration tubing and two holders for fixing the tubing to one of the sensor rods. 
6.2769.120	EMPTY sensor for level sensor for 10 L canisters	This sensor, mounted in the 6.1626.000 bottle attachment, results in the level sensor EMPTY for 10 L canisters. Included are two sensor rods, one aspiration tubing and two holders for fixing the tubing to one of the sensor rods. 