



# 891 Professional Analog Out

2.891.0010

The 891 Professional Analog Out is an additional module for intelligent IC instruments. It allows to transfer analog signals of different chromatographic data to external data systems.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

## Scope of delivery 2.891.0010

| Qt.   | Order no.  | Description                 |
|-------|------------|-----------------------------|
| 1 PCS | 1.891.0010 | 891 Professional Analog Out |

The 891 Professional Analog Out is an additional module for intelligent IC instruments. It allows to transfer analog signals of different chromatographic data to external data systems.



1 PCS

6.2128.200

Cable 891-Analog

Analog cable to connect to external AD converters.



---

1 PCS

6.2156.070

Cable Extension Module - Professional IC, 1 m

Cable connecting the extension module with the Professional IC instrument

