



NIRS XDS Monochromator

2.921.0010

The result of 50 years experience in the area of spectroscopy was incorporated in the development of this NIRSystems spectrometer. The unique technology improves the signal-noise ratio at high resolution and enhances accuracy and precision. The patented monochromator and the modular design of the measuring module system make NIRSystems instruments indispensable for the spectral measurement of every imaginable sample type in the visible and near-infrared wavelength range (400 - 2,500 nm). The measuring modules can be replaced without difficulty, even during operation. This means that the user profits from the great flexibility in all areas of application.

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 2.921.0010

Qt.	Order no.	Description
-----	-----------	-------------

The result of 50 years experience in the area of spectroscopy was incorporated in the development of this NIRSystems spectrometer. The unique technology improves the signal-noise ratio at high resolution and enhances accuracy and precision. The patented monochromator and the modular design of the measuring module system make NIRSystems instruments indispensable for the spectral measurement of every imaginable sample type in the visible and near-infrared wavelength range (400 - 2,500 nm). The measuring modules can be replaced without difficulty, even during operation. This means that the user profits from the great flexibility in all areas of application.



Accessory kit for NIRS XDS laboratory monochromators

