



Metrosep Chel PCC 1 VHC/4.0

6.01010.350

Preconcentration column for polyvalent cations. Monovalent cations are removed by matrix elimination prior to injection.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Scope of delivery 6.01010.350

Qt.	Order no.	Description
2 PCS	6.2744.060	Threaded stopper
		For UNF 10/32. Stopper for IC, e.g. for the sealing of columns.

