




Titration head 2xNS14 (P&P)

6.01403.040

Titration head with 2 NS14 openings for use in Pick&Place modules of the OMNIS Sample Robot Pick&Place with 3 retracted dosing tips. Prepared for being equipped with rinsing nozzles. One aspiration tube is already retracted into the titration head.

Below, the accessories are grouped into Scope of delivery and Optional accessories. Please keep this printout at hand for ordering replacement material. These lists may be subject to change.

Scope of delivery 6.01403.040

Qt.	Order no.	Description
2 PCS	6.1236.020	Sleeve with SGJ 14/12 mm
		Sleeve with SGJ 14/12 mm and O-ring.



1 PCS

6.1446.000

SGJ stopper / B-14/(15)



4 PCS

6.1808.000

Connection bushing M6



1 PCS

6.1808.230

Adapter M6 / M8 inner

Adapter with M6 and M8 inner threads



3 PCS

6.2726.090

Antidiffusion valve

For antidiffusion buret tips 6.1541.XXX, 6.1543.130, 6.1543.140,
6.1543.150, 6.1543.160, 6.1543.200



Optional accessories

Order no.	Description	
6.1829.050	Replacement tubings for retracted tubings in titration heads and robotic arms Three replacement tubings with tapering tip.	
