



Titration head 3xNS14 / 4x6.4 mm for 200 mL disposable beaker (P&P)

6.01403.060

Titration head with 3 NS14 and 4x 6.4 mm openings for using 200 mL PP disposable beakers (6.1459.310) or 150 mL beaker glasses (high shape, without spout) in Pick&Place modules of the OMNIS Sample Robot Pick&Place with two retracted dosing tips. Prepared for being equipped with rinsing and aspiration tubing.

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.


Optional accessories

Order no.	Description
-----------	-------------

6.1829.050

Replacement tubings for retracted tubings in titration heads and robotic arms

Three replacement tubings with tapering tip.



6.1459.310

Sample beaker, disposable, plastic (PP), 200 mL, 1000 pieces

Sample beaker, disposable, plastic (polypropylene, PP), 200 mL, suitable for 6.02041.020 OMNIS Software sample rack, 6.2041.370 und 6.2041.830 sample racks, 1000 pieces.

