



5 x SGJ 14 titration vessel lid OMNIS

6.01405.000

Titration vessel lid with five SGJ 14 openings for 601406XX0 titration vessels with screw cap.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Optional accessories

Order no.	Description
6.1440.010	Gas inlet and overflow tube with valve For bubbling or covering of the sample by inert gas.



6.2727.010

Intermediate SGJ sleeve

Used with reference and auxiliary electrodes in measuring and titration vessel with 6.1414.010 lid

