



OMNIS titration vessel lid for volumetric Karl Fischer titration

6.01405.010

Titration vessel lid for volumetric Karl Fischer titration with an OMNIS titrator.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Scope of delivery 6.01405.010

Qt.	Order no.	Description
1 PCS	6.1454.130	Form seal made of FPM for 61414xxx
		for KF and VA instruments



Optional accessories

Order no.	Description	
6.02709.030	Screw nipple Plastic screw nipple, suitable for 6014050X0 titration vessel lids.	
6.02709.020	Septum stopper Plastic septum stopper, suitable for 6014050X0 titration vessel lids.	
6.02709.040	Stopper Stopper for titration vessel lids 6014050X0.	